



70375

PIPES Buffer extrapure, 99%

Part B

CAS : 5625-37-6

Molecular Formula : C₈H₁₈N₂O₆S₂

Molecular Weight : 302.36

Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Solubility (Turbidity) 5% solution in 10% NaOH	Clear
Solubility (Colour) 5% solution in 10% NaOH	Colourless
Assay (T)	min. 99%
pH (0.1% aq. solution)	2.5 - 3.0
Loss on drying	max. 0.5%
Sulphated Ash	max. 0.1%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.001%
Absorbance (A) of 10% in 10% NaOH in a 1cm cell @250 nm	max. 0.2

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29335990 (GST 18%)
25 Gms	29335990 (GST 18%)
500 Gms	29335990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
