



70427

Tetradecylammonium Bromide (TDAB) extrapure AR, 99%

Part A

CAS : 14937-42-9

Molecular Formula : C₄₀H₈₄BrN

Molecular Weight : 659.01

Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Solubility (Turbidity) 5% solution in Methanol	Clear
Solubility (Colour) 5% solution in methanol	Colourless
Assay (AT)	min. 99%
Melting Point	87 - 90°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29239000 (GST 18%)
25 Gms	29239000 (GST 18%)
5 Gms	29239000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com