



70514

## Lauric Acid pure, C12-99%

Part A

CAS : 143-07-7

Molecular Formula : C<sub>12</sub>H<sub>24</sub>O<sub>2</sub>

Molecular Weight : 200.32

### Specifications

Appearance (Colour)	White
Appearance (Form)	Flakes or granules
Assay (GC)	min. 99%
Iodine value	max. 0.30
Saponification value	279-283
Acid value	278-282

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2.5 Kg	29159010 (GST 18%)
500 Gms	29159010 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

### Available Packages

500 Gms  
5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com