



70851

## Sodium Molybdate Dihydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 10102-40-6

Molecular Formula : Na<sub>2</sub>MoO<sub>4</sub>.2H<sub>2</sub>O

Molecular Weight : 241.95

### Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (T)	min. 99%
pH (5% aq. solution)	7.0 - 10.5
Water Insoluble matter	max. 0.005%
Chloride (Cl)	max. 0.001%
Sulphate (SO <sub>4</sub> )	max. 0.015%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.0005%
Phosphate (PO <sub>4</sub> )	max. 0.0005%
Nitrate (NO <sub>3</sub> )	max. 0.005%
Ammonium (NH <sub>4</sub> )	max. 0.001%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28417020 (GST 18%)
250 Gms	28417020 (GST 18%)
500 Gms	28417020 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
250 Gms	Plastic Bottle
500 Gms	Plastic Bottle

## Available Packages

100 Gms

250 Gms

500 Gms

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---