



70874

Triphenylphosphonium Dibromide pure, 98%

Part B

CAS : 1034-39-5

Molecular Formula : C₁₈H₁₅Br₂P

Molecular Weight : 422.09

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
HSN Code	
100 Gms	29319090 (GST 18%)
25 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)

Available Packages

5 Gms

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com