



70940

## Anisole extrapure, 99%

Part B

CAS : 100-66-3

Molecular Formula : C<sub>7</sub>H<sub>8</sub>O

Molecular Weight : 108.14

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.993-0.995
Refractive Index (20°C)	1.516-1.518
Boiling Range	152-156°C
Phenol (GC)	max. 0.2%
Water (KF)	max. 0.3%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2222 , IMCO Class No.:3 , Packing Group:III
HSN Code	
500 ml	29093019 (GST 18%)

### Available Packages

500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com