



71099

Crucible Porcelain, 30ml, Wide with Lid

Part L

Specifications

Appearance (Colour)	White to off White
Appearance (Form)	Crucible
Material of Construction (MOC)	Porcelain
Crucible Dia	Top OD 50mm \pm 2
Crucible Outer Height	31mm \pm 3
Crucible Wall Thickness	1.8mm \pm 0.9
Crucible Weight	30gm \pm 5
Max Usable Temperature	1150°C \pm 50
Flextural Strength @20°C	min.80 Mpa (815Kg/cm ²)
Lid Dia	ID 53mm \pm 5

General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
4 Nos	69091100 (GST 18%)

Available Packages

4 Nos

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com