



71548

Silica Gel Orange (Self Indicating), 2-5mm

Part A

CAS : 112926-00-8

Molecular Formula : SiO₂

Specifications

Appearance (Colour)	Orange
Appearance (Form)	Beads/Crystalline compound
Loss on drying	max. 5%
Mesh size	5 - 8
Water absorption capacity	min. 25%

Other Information

Notes

This type of Silica gel, is impregnated with an "organic indicator" which is rated as completely safe. When free from moisture, the particles are dark orange. As the particles take up moisture, they turn to Orange-Yellow and then to Green and this indicates saturation and now requires replacement or regeneration.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	28112200 (GST 18%)
500 Gms	28112200 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com