



71772

## Molybdic Acid extrapure AR, ACS, ExiPlus, Multi-Compendial, 85%

Part A

CAS : 7782-91-4

### Specifications

ExiPlus	meets compendial specs of USP
Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in dilute Ammonia	Clear
Solubility (Colour) 10% solution in dilute Ammonia	Colourless
Assay (T)	min. 85%
Insoluble in dilute Ammonium Hydroxide	max. 0.01%
Chloride (Cl)	max. 0.002%
Arsenate, Phosphate and Silicate (PO4)	max. 0.001%
Phosphate (PO4)	max. 0.0005%
Sulphate (SO4)	max. 0.02%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	28257020 (GST 18%)
100 Gms	28257020 (GST 18%)
250 Gms	28257020 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com