



71801

**Decane Sulphonic Acid Sodium Salt
Monohydrate for HPLC, 99%**

Part A

CAS : 13419-61-9

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (Ion exchange)	min. 99%
pH (10% aq. solution)	6.0 - 7.0
Loss on drying at 105°C	6.0 - 7.0 %
Absorbance of 0.005M in water,1cm cell vs HPLC grade water	
@200nm	max.0.15
@220nm	max. 0.05
Absorbance of 0.1M in water,1cm cell vs HPLC grade water	
@254nm	max.0.05

General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29041090 (GST 18%)
25 Gms	29041090 (GST 18%)
250 Gms	29041090 (GST 18%)

Available Packages

25 Gms
100 Gms
250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com