



72033

## Copper Metal powder extrapure, 99.5%, -200 mesh

Part A

CAS : 7440-50-8

Molecular Formula : Cu

Atomic Weight : 63.55

### Specifications

Appearance (Colour)	Fine reddish brown
Appearance (Form)	Metallic powder
Solubility (Turbidity) 10% solution in 5N HNO <sub>3</sub>	Clear
Solubility (Colour) 10% solution in 5N HNO <sub>3</sub>	Blue
Assay	min. 99.5%
Particle Size	~200 mesh BSS

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3089 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
500 Gms	74061000 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com