



72217

Bismuth Metal Powder extrapure, 99%, -325 mesh

Part A

CAS : 7440-69-9

Molecular Formula : Bi

Atomic Weight : 208.98

Specifications

Appearance (Colour)	Greyish black
Appearance (Form)	Powder
Solubility (Turbidity) 5% solution in 5N HNO3	Clear
Assay (EDTA Titration)	min. 99%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	81061020 (GST 18%)
500 Gms	81061020 (GST 18%)

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com