



72367

**4-Nitro-L-Phenylalanine Monohydrate extrapure,  
99%**

Part B

CAS : 207591-86-4

Molecular Formula : C<sub>9</sub>H<sub>10</sub>N<sub>2</sub>O<sub>4</sub>.H<sub>2</sub>O

Molecular Weight : 228.20

**General Information**

Storage	2 to 8°C (Refrigerate)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29420090 (GST 18%)

**Available Packages**

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com