



72918

Bismuth (III) Nitrate Pentahydrate extrapure, 98%

Part A

CAS : 10035-06-0

Molecular Formula : Bi(NO3)3.5H2O

Molecular Weight : 485.07

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% solution in 5N HNO3	Clear
Solubility (Colour) 10% solution in 5N HNO3	Colourless
Assay (EDTA Titration)	min. 98%
Chloride (Cl)	max.0.005%
Sulphate (SO4)	max.0.02%
Sodium (Na)	max. 0.02%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1477 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
100 Gms	28342990 (GST 18%)
500 Gms	28342990 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com