



73017

Immersion Oil Synthetic for microscopy

Part A

CAS : 8000-27-9

Specifications

Appearance (Colour)	Colourless to pale yellow
Appearance (Clarity)	Clear
Appearance (Form)	Viscous liquid
Refractive Index (20°C)	1.515 - 1.516
Density (g/ml) @ 20°C	1.10 - 1.30 g/ml

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
125 ml	39100090 (GST 18%)
30 ml	39100090 (GST 18%)
5 Ltr	39100090 (GST 18%)
500 ml	39100090 (GST 18%)

Available Packages

30 ml

125 ml

500 ml

5 Ltr

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com