



73081

Silicon Metal powder, 98.5%, -200 mesh

Part B

CAS : 7440-21-3

Atomic Formula : Si

Atomic Weight : 28.09

Specifications

Appearance (Colour)	Greyish Black
Appearance (Form)	Powder
Assay	min. 98.5%
Mesh size	-200

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1346 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
100 Gms	28046900 (GST 18%)
500 Gms	28046900 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com