



73215

## Acetic Acid Dried, 99.9%, water 0.05%

Part A

CAS : 64-19-7

Molecular Formula : C<sub>2</sub>H<sub>4</sub>O<sub>2</sub>

Molecular Weight : 60.05

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max.10
Assay (GC)	min. 99.9%
Density (g/ml) @ 20°C	1.048-1.049
Freezing Point	~16°C
Boiling Point	118°C
Non Volatile Matter	max. 0.001%
Chloride (Cl)	max. 0.0001%
Sulphate (SO <sub>4</sub> )	max. 0.0001%
Iron (Fe)	max. 0.00002%
Heavy Metals (Pb)	max. 0.000002%
Substances reducing Dichromate (O)	max. 0.003%
Formate (HCOO)	max. 0.01%
Water (KF)	max. 0.05%

### Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2789 , IMCO Class No.:8(3) , Packing Group:II
HSN Code	
1000 ml	29152100 (GST 18%)

500 ml

29152100 (GST 18%)

Type of Packing

1000 ml

Glass Bottle

500 ml

Glass Bottle

### Available Packages

500 ml

1000 ml

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---