



73229

o-Phenylenediamine Dihydrochloride (OPD.2HCl) extrapure AR, 99%

Part B

CAS : 615-28-1

Molecular Formula : C₆H₈N₂.2HCl

Molecular Weight : 181.06

Specifications

Appearance (Colour)	White to pinkish to tan to grey
Appearance (Form)	Crystalline powder
Assay	min. 99%
Loss on drying	max. 0.2%

Other Information

Application	Commonly used ELISA substrate.
-------------	--------------------------------

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2811 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
25 Gms	29215990 (GST 18%)
5 Gms	29215990 (GST 18%)

Available Packages

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com