



73342

Nickel Aluminium Alloy powder, 50% Ni, 50% Al

Part A

CAS : 12635-27-7

Specifications

Appearance (Colour)	Grey - black
Appearance (Form)	Powder
Nickel (Ni) content	~50%
Aluminium (Al) content	~50%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3089 , IMCO Class No.:4.1 , Packing Group:II
HSN Code	
100 Gms	75022030 (GST 18%)
500 Gms	75022030 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com