



73694

2-Pyrrolidone pure, 99%

Part B

CAS : 616-45-5

Molecular Formula : C₄H₇NO

Molecular Weight : 85.11

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Brown
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 25°C	1.109-1.113gm
Refractive Index (20°C)	1.486-1.487
Melting Point	23-25°C
Water	max. 0.3%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 ml	29337990 (GST 18%)
2500 ml	29337990 (GST 18%)
500 ml	29337990 (GST 18%)
Type of Packing	
250 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

250 ml

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com