



73698

Fluorescein Isothiocyanate Isomer I (FITC) extrapure, 95%

Part B

CAS : 3326-32-7

Molecular Formula : C₂₁H₁₁NO₅S

Molecular Weight : 389.39

Specifications

Appearance (Colour)	Yellow to Orange to Dark Orange Red
Appearance (Form)	Powder
Absorption maxima	488 - 495 nm (0.1 N NaOH)
Assay (UV)	min. 95%

Other Information

Application

Applications:
- Used for conjugation to all of our antibodies and purified proteins, with the exception of streptavidin.

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	32049000 (GST 18%)
100 Mg	32049000 (GST 18%)
25 Mg	32049000 (GST 18%)
250 Mg	32049000 (GST 18%)
Type of Packing	
100 Mg	Glass Bottle
25 Mg	Glass Bottle
250 Mg	Glass Bottle

Available Packages

25 Mg

100 Mg

250 Mg

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
