



73932

## Zinc Metal Dust extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%, -325 mesh

Part A

CAS : 7440-66-6

Molecular Formula : Zn

Atomic Weight : 65.38

### Specifications

ExiPlus	Meets compendial specs of Ph.Eur
Appearance (Colour)	Grey
Appearance (Form)	Fine powder free from aggregates
Solubility (Turbidity) 10% solution 5N HCl	Clear
Solubility (Colour) 10% solution 5N HCl	Colourless
Assay (EDTA Titration)	min. 99.5%
Iron (Fe)	max. 0.005 %
Copper (Cu)	max. 0.0005 %
Lead (Pb)	max. 0.01%
Fatty matter	Not detected
Substances reducing Permanganate (O)	max. 0.0015%
Acid insoluble matter	max. 0.01%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1436 , IMCO Class No.:4.3(4.2) , Packing Group:II
HSN Code	
2.5 Kg	79031000 (GST 18%)
500 Gms	79031000 (GST 18%)

### Available Packages

500 Gms

2.5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com