



73961

Eschka's Mixture extrapure AR, 38%

Part A

CAS : 8007-09-8

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline Powder
Solubility (Turbidity) 10% solution 5N HCl	Clear
Magnesium (Mg) content	~ 38%
Loss on drying	max. 1%
Sulphate & sulphite (SO ₄)	max 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	38229090 (GST 12%)
250 Gms	38229090 (GST 12%)
Type of Packing	
1 Kg	Plastic Bottle

Available Packages

250 Gms

1 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com