



74118

Thioglycolic Acid ExiPlus, Multi-Compendial, 80%

Part A

CAS : 68-11-1

Molecular Formula : C₂H₄O₂S

Molecular Weight : 92.12

Specifications

ExiPlus	meets compendial specs of USP
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (by Neutralization)	~80%
Density (g/ml) @ 20°C	1.265-1.275
Refractive Index (20°C)	1.471-1.472
Identification Test	Passes test
Sensitivity for Iron (Fe)	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1940 , IMCO Class No.:8 , Packing Group:II
HSN Code	
2500 ml	29309099 (GST 18%)
500 ml	29309099 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com