



74193

## Paraffin Wax White Soft pure

Part A

CAS : 8002-74-2

### Specifications

Appearance (Colour)	White
Appearance (Form)	Translucent soft mass
Melting range	38 - 56°C
Acidity	Passes test
Sulphated Ash	max. 0.1%

### Other Information

Application

Paraffin Wax is a mixture of hydrocarbons containing between 20-40 carbon atoms that are derived out of petroleum extractions. It is a white, odourless, tasteless waxy solid which is used in lubrication and histology applications.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	27129090 (GST 18%)
500 Gms	27129090 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com