



74212

Tetrabutylammonium Hydroxide 40% in Methanol extrapure

Part A

CAS : 2052-49-5

Molecular Formula : C₁₆H₃₇NO

Molecular Weight : 259.48

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (T)	~40%
Halides (Cl)	max 0.04%
Sulphate (SO ₄)	max 0.3%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3267 , IMCO Class No.:8 , Packing Group:II
HSN Code	
100 ml	29239000 (GST 18%)
250 ml	29239000 (GST 18%)
500 ml	29239000 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

100 ml

250 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com