



74380

Perovskite Quantum Dots (pQD 450), Fluorescence 450nm (Hydrophobic)

Part C

Specifications

Appearance (Colour)	Yellowish green
Appearance (Form)	Crystals
Emission wavelength	450nm

Other Information

Notes

Readily soluble in toluene, chloroform and similar solvents. Not soluble in alcohols, ethers and polar solvents.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	36 Months
HSN Code	
50 Mg	28419000 (GST 18%)

Available Packages

50 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com