



74384

Methyl Methanesulphonate (MMS) extrapure, 99%

Part B

CAS : 66-27-3

Molecular Formula : C₂H₆O₃S

Molecular Weight : 110.113

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.298 - 1.300 g/ml
Refractive Index (20°C)	1.413 - 1.414

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2810 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 ml	98020000 (GST 18%)
25 ml	98020000 (GST 18%)

Available Packages

25 ml

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com