



74495

Uranyl Nitrate Hexahydrate extrapure, 98%

Part B

CAS : 13520-83-7

Molecular Formula : $UO_2(NO_3)_2 \cdot 6H_2O$

Molecular Weight : 502.13

Specifications

Appearance (Colour)	Yellow
Appearance (Form)	Compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Yellow
Assay	min. 98%
Melting Point	~ 60°C

Other Information

Application

Uranium nitrate is used as a negative stain for viruses in electron microscopy; in tissue samples it stabilizes nucleic acids and cell membranes.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:1477 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
1 Gms	28443029 (GST 18%)
5 Gms	28443029 (GST 18%)
Type of Packing	
1 Gms	Plastic Bottle
5 Gms	Glass Bottle

Available Packages

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com