



74710

## 1-Docosanol (Behenyl Alcohol) extrapure, 98%

Part B

CAS : 661-19-8

Molecular Formula : C<sub>22</sub>H<sub>46</sub>O

Molecular Weight : 326.61

### Specifications

Appearance (Colour)	White
Appearance (Form)	Powder or Pellets/Flakes
Assay (GC)	min. 98%
Melting Point	65 - 72° C

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29051990 (GST 18%)
500 Gms	29051990 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

### Available Packages

25 Gms

---

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com