



74726

MacConkey Broth (MacConkey Broth w/ Neutral Red)

Part D

Other Information

Applications

For the selective isolation and cultivation of coliforms in milk and water.

Composition

Ingredients	gm/lt.
Peptone	20.00
Lactose	10.00
Sodium chloride	5.00
Bile salts	5.00
Neutral Red	0.075

Directions

1. Add 40.00 gm powder to 1.0 liter of distilled/purified water and mix thoroughly.
2. Gently heat and bring to boiling.
3. Dispense into tubes containing inverted Durham's tubes.
4. Autoclave at 15 psi pressure at 121°C for 15 minutes.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	38210000 (GST 18%)
500 Gms	38210000 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com