



74926

D(+)-Arabitol (D-Arabinitol) extrapure 99%

Part B

CAS : 488-82-4

Molecular Formula : C₅H₁₂O₅

Molecular Weight : 152.15

Specifications

Appearance (Colour)	Off-white to faint beige
Appearance (Form)	Powder
Solubility (Turbidity) 5% aq. solution	Clear
Assay	min. 99%
Melting Point	101 - 104°C
Specific rotation [?]	+10° to +12° (C=5 , 8% Borax soln.)

Other Information

Application	Substrate for D-Arabitol Dehydrogenase
-------------	---

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29400000 (GST 18%)
25 Gms	29400000 (GST 18%)
5 Gms	29400000 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

1 Gms

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com