



75122

## Xylene Cyanol FF extrapure AR, ExiPlus, Multi-Compendial

Part A

CAS : 2650-17-1

Molecular Formula : C<sub>25</sub>H<sub>27</sub>N<sub>2</sub>NaO<sub>6</sub>S<sub>2</sub>

Molecular Weight : 538.6

### Specifications

ExiPlus	meets compendial specs of USP
Appearance (Colour)	Green
Appearance (Form)	Crystalline powder
Loss on drying	max. 8%
Absorbance (A) of 1% aq. solution in a 1cm cell v/s H <sub>2</sub> O	
@614nm	1000 - 1400
Molar Extinction coefficient (?) in Water at 613 - 616 nm	min. 65000 (On dried basis)

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	32041248 (GST 18%)
25 Gms	32041248 (GST 18%)
5 Gms	32041248 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

### Available Packages

5 Gms

25 Gms

100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com