



75241

2-Thiobarbituric Acid (TBA) extrapure AR, 99%

Part A

CAS : 504-17-6

Molecular Formula : C4H4N2O2S

Molecular Weight : 144.15

Specifications

Appearance (Colour)	White to off-white to cream to pale yellow
Appearance (Form)	Powder
Solubility (Turbidity) 5% solution in 1N NaOH	Clear to slight hazy
Solubility (Colour) 5% solution in 1N NaOH	Almost colourless to pale yellow
Assay (T) (on dry basis)	min. 99.0%
Water (KF)	max. 1.0%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3335 , IMCO Class No.:9 , Packing Group:III
HSN Code	
100 Gms	29335200 (GST 18%)
25 Gms	29335200 (GST 18%)

Available Packages

25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com