



75311

1-H-Tetrazole sublimed extrapure, 99%

Part B

CAS : 288-94-8

Molecular Formula : CH₂N₄

Molecular Weight : 70.06

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 2.5% solution in Acetonitrile	Clear
Solubility (Colour) 2.5% solution in Acetonitrile	Colourless
Assay	min 99%
Melting Point	155 - 157°C

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:1993 , IMCO Class No.:3 , Packing Group:II
HSN Code	
1 Gms	29331999 (GST 18%)
25 Gms	29331999 (GST 18%)
5 Gms	29331999 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle

Available Packages

1 Gms

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com