



75531

BCIP Red/NBT Solution A (BCIP Red Solution)

Part B

Specifications

Appearance (Clarity)	Clear
Appearance (Form)	Colourless to pale brown liquid
Composition	1% solution of 5-bromo-4-chloro-3-indolylphosphate-p-Toluidine salt (BCIP) in N,N-dimethylformamide (DMF)
Suitability as chromogenic substrate	Passes

Other Information

Application

Applications:
It is a chromogenic substrate for enzyme alkaline phosphatase, producing a blue coloured precipitate at the site of enzymatic activity. It has applications in immunohistochemistry and immunocytochemistry.

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	12 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	38229090 (GST 12%)

Available Packages

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com