



75543

DABCYL-N-Succinimidyl Ester extrapure, 98%

Part B

CAS : 146998-31-4

Molecular Formula : C19H18N4O4

Molecular Weight : 366.37

Specifications

Appearance (Colour)	Orange red
Appearance (Form)	Powder
Assay	min. 98%

Other Information

Application

A useful reagent (similar to MANCYL-SE) for labeling proteins or peptides through their amino-groups by forming stable peptide bonds. The succinimidyl ester is reactive with terminal amines or lysines of peptides and other nucleophiles for fluorescent studies of proteins. Its fluorescent properties (has a characteristic broad and intense visible absorption but has no fluorescence) make it an ideal long wavelength quencher.

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
HSN Code	
250 Mg	29337990 (GST 18%)
50 Mg	29337990 (GST 18%)
Type of Packing	
250 Mg	Plastic Vial
50 Mg	Plastic Vial

Available Packages

50 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.  
608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com