



76043

Sterility Test Kit (Pre-Weighed Media)

Part D

Specifications

Content

Kit contains..
(46990) Fluid Thioglycollate
Medium,Harmonized 29.75 gm/ltr
(10935) Soyabean Casein Digest
Medium,(SCDM),Tryptone Soya
Broth,(TSB),Harmonized. 30.00 gm/ltr.

Other Information

Applications

Used for sterility testing of products as per Pharmacoepl standards.

Composition

Fluid Thioglycollate Medium, Harmonized	29.75
Soyabean Casein Digest Medium,(SCDM),Tryptone Soya Broth,(TSB),Harmonized	30.00 gm/lt.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com