



76072

Creatine Monohydrate extrapure, 99%

Part B

CAS : 6020-87-7

Molecular Formula : C₄H₉N₃O₂.H₂O

Molecular Weight : 149.15

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay (NT)	min. 99%
Loss on drying	11.8 - 12.4%
Chloride (Cl)	max. 0.01%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%
Creatinine	max. 0.2%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	29224990 (GST 18%)
100 Gms	29224990 (GST 18%)
25 Gms	29224990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

Available Packages

25 Gms
100 Gms
1 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com