



76171

Agarose TBE mix 1.0%

Part B

Other Information

Notes

Pre-mixed blend of Agarose and TBE
at 1x concentration.

General Information

Storage

25 to 40°C (Room Temperature)

IMDG Identification

Not Regulated for Transport (Non-Haz)

HSN Code

100 Gms

38229090 (GST 12%)

Type of Packing

100 Gms

Plastic Bottle

Available Packages

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com