



76207

Oxidase Reagent (Gordon-MacLeod Reagent)

Part B

Specifications

Appearance (Colour)	Black
Appearance (Form)	Liquid
Suitability for Microscopy	Passes test

Other Information

Application	The reagent is used for detection of oxidase activity of a bacterial culture.
-------------	---

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	38229090 (GST 12%)

Available Packages

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com