



76265

1-Tetralone pure, 98%

Part B

CAS : 529-34-0

Molecular Formula : C₁₀H₁₀O

Molecular Weight : 146.19

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Brownish yellow
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	1.096-1.098
Refractive Index (20°C)	1.568-1.569
Water	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	29146990 (GST 18%)
500 ml	29146990 (GST 18%)
Type of Packing	
100 ml	Glass Bottle

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com