



76333

## 1,1,2,2-Tetrachloroethane extrapure AR, 99%

Part B

CAS : 79-34-5

Molecular Formula : C<sub>2</sub>H<sub>2</sub>Cl<sub>4</sub>

Molecular Weight : 167.85

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 20
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.596-1.597
Refractive Index (20°C)	1.493-1.494
Boiling Range	144-146°C
Non Volatile Matter	max. 0.005%
Acidity (HCl)	max. 0.002%
Substance Discoloured by H <sub>2</sub> SO <sub>4</sub>	Passes test
Chloride (Cl)	max. 0.0002%
Iron (Fe)	max. 0.00005%
Heavy Metals (Pb)	max. 0.0001%
Free Chlorine (Cl)	max. 0.0005%
Water (KF)	max. 0.05%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1702 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
1000 ml	29031910 (GST 18%)
500 ml	29031910 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle
500 ml	Glass Bottle

## Available Packages

500 ml

1000 ml

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---