



76512

Eriochrome Cyanine R

Part A

CAS : 3564-18-9

Molecular Formula : C₂₃H₁₅Na₃O₉S

Molecular Weight : 536.40

Specifications

Appearance (Colour)	Brick red
Appearance (Form)	Powder
Absorbance (A) of 1% aq. solution (pH 7.0) in a 1cm cell	
@438nm	130 - 200
Loss on drying	max.10%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	32041990 (GST 18%)
25 Gms	32041990 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

10 Gms

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com