



76601

Starch Soluble extrapure AR, ACS, ExiPlus, Multi-Compendial

Part A

CAS : 9005-84-9

Specifications

ExiPlus	meets compendial specs of USP
Appearance (Colour)	White
Appearance (Form)	Powder
Solubility (Turbidity) 2% aq. hot solution	Slight opalescent
Solubility (Colour) 2% aq. hot solution	Colourless
pH (2% aq. solution)	5.0 - 7.0
Loss on drying	max. 15%
Chloride (Cl)	max. 0.005%
Sulphated Ash	max. 0.4%
Acidity	max. 0.8ml N%
Dextrin	Passes test
Reducing sugars (C ₆ H ₁₂ O ₆)	max. 0.2%
Sensitivity to iodine	Passes test
Solution quality	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	35051090 (GST 18%)
500 Gms	35051090 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com