



76610

Silver Oxide extrapure AR, 99%

Part A

CAS : 20667-12-3

Molecular Formula : Ag₂O

Molecular Weight : 231.74

Specifications

Appearance (Colour)	Greyish Black
Appearance (Form)	Powder
Solubility (Turbidity) 1% solution in 5N HNO ₃	Clear to Slight opalascnt
Assay (AT)	min. 99%
Alkalinity	max. 5ml N%
Nitrate (NO ₃)	max. 0.05%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1479 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
10 Gms	28432900 (GST 18%)
25 Gms	28432900 (GST 18%)

Available Packages

10 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com