



77031

Biuret Reagent Solution ExiPlus, Multi-Compendial

Part B

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Blue
Appearance (Form)	Liquid
Suitability as indicator for Protein	Passes test

Other Information

Application	Used for protein determination.
-------------	---------------------------------

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3266 , IMCO Class No.:8, Packing Group:III
HSN Code	
250 ml	38229090 (GST 12%)

Available Packages

250 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com