



78084

Popular range Protein Marker: 14-66 kDa

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Blue
Appearance (Clarity)	Clear
SDS-PAGE Analysis	Electrophoresis of the marker should show six crisp bands of almost equal intensity and be visualised by Coomassie Blue staining or SpryBlue staining.

Other Information

Description	BioLit Popular range Protein Marker covers the typical range of protein sizes that a researcher or student requires for routine experiments. This is a ready-to-use protein marker. The markers have been carefully prepared under sterile conditions and stored in appropriate additives to ensure long term stability.
Application	SDS PAGE Analysis & Western Blotting

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	18 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
0.5 ml	38229090 (GST 12%)
1 ml	38229090 (GST 12%)

Available Packages

0.5 ml
1 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com