



78339

Silica C18 (17%) Bonded monomeric extrapure

Part B

Specifications	
Appearance (Colour)	White
Appearance (Form)	Powder
Particle Size	40 - 63 µm
Pore size	60 Å

Other Information	
Application	This versatile Silica C18 is the best choice for your applications in reversed-phase chromatography.

General Information	
Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
HSN Code	
25 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)
Type of Packing	
25 Gms	Plastic Bottle
5 Gms	Plastic Bottle

Available Packages	
5 Gms	
25 Gms	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com