



78523

L-Homophenylalanine extrapure, 97%

Part B

CAS : 943-73-7

Molecular Formula : C₁₀H₁₃NO₂

Molecular Weight : 179.22

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min. 97%
Specific rotation [?]	+ 45° to + 46° (C=2, 3M HCl)

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29224990 (GST 18%)
25 Gms	29224990 (GST 18%)
5 Gms	29224990 (GST 18%)
500 Mg	29224990 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

500 Mg

1 Gms

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com